La almohada-soporte ortopédica

Asi se duerme con The Pillow

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The orthopedic pillow

How to Sleep with The Pillow

The Pillow is an orthopedic pillow designed by Dr. John Downing (U.S.A.) to support healthy sleep posture. It's important to use it the right way. So please take a moment to read these instructions carefully. The Pillow consists of a neck roll (1) at the front, two side supports (2), and a thinner, flexible center section (3).


Sleeping on Your Back

When you lie on your back, your head rests on the flexible center (3). The roll under your shoulders creates a smooth transition from shoulder to neck to head. This avoids pinching or pressure in your neck. In this position, your body can fully relax. Common pressure points that other pillows cause are gently relieved.

Sleeping on Your Side

Sleeping on your side or shifting from your side to your back feels completely natural with The Pillow. It stays in place as you move, gently guiding your head onto the slightly raised side sections that provide just the right support.

Your body may need a little time to adjust to this healthier way of sleeping. That's completely normal. Be patient and give yourself a few nights. Before you know it, you'll wonder how you ever slept without The Pillow. Some people take a bit longer to get used to it, while others are surprised by how comfortable it feels right from the start. To keep The Pillow in top condition, avoid placing it in direct sunlight or bright daylight.

The orthopedic pillow for headaches, poor sleep, and neck or back pain.

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